



# Antibak POWDER

High Level  
Disinfectant Powder

# Antibak POWDER



A multipurpose High Level Disinfectant and cleaner perfect for both large and small scale cleaning & decontamination. Easy to use sachets guarantee the correct dilution every time - without the need for any special PPE!

## Overview

### Meets International Standards

Tested against over 30 pathogens and shown to be effective against viruses, spores, fungi, bacteria and mycobacteria

### Easy to use sachet

No mess or fuss – correct dilution every time!

### Versatile

Compatible with all surfaces, when used as instructed

### Non-hazardous solution

When diluted, AntiBak Powder requires no hazard labelling - safer for you and the environment

### Long lasting solution

Solution remains active for up to 7 days

### Safe and easy disposal

Can be disposed of to the drain as solution degrades to harmless salts

## How to use

50g sachets makes 5 litres of solution. For use in open top buckets.

- To use at a 1% solution, add the entire contents of one 50g sachet to 5 litres of water.
- Always use entire sachet contents - never store opened sachets.

### Step 1



Open sachet.

### Step 2



Pour powder onto the surface of the water.

### Step 3



Stir for up to two minutes before use.

### Step 4



Apply the solution using a mop & bucket, cloths, wipes or other cleaning equipment.

### Step 5



Rinse surfaces after application, if applicable. Rinse all cleaning equipment in clean water after use.

## Where to use

AntiBak Powder is designed for the deep cleaning and High Level Disinfection of areas and items such as:

- Isolation rooms
- Laboratories
- Mortuaries
- Sluice rooms
- Walls & floors
- Blood spills
- Bodily fluid spills
- Steel gurneys
- Door handles
- Beds & mattresses
- Commodes
- HVAC
- Toilets, sinks & taps
- Dressing trolleys
- Macerators

## Compatibility

Tested for compatibility against materials including:

- Aluminium
- Steel
- Rubber
- Pyrex
- Acrylic
- PVC
- Mattress fabric & covers

Not recommended for use on surfaces containing tin, copper, brass or zinc, or on polished timber.

## Pack sizes

Description	Case Size	Order Code
50g sachet	300 per case	ABP-S-50

ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE. ALWAYS FOLLOW CLEANING PROCEDURES AND GUIDELINES.

## Test results

All results are from use at a 1% dilution.

### Viruses

Pathogen	Test	Pass time
Adenovirus	EN 14476	15 min.
Bovine Herpes Virus 1	Independent <sup>1</sup>	5 min.
Bovine Viral Diarrhoea Virus	Independent <sup>2</sup>	5 min.
Feline Calici Virus	Independent <sup>1</sup>	5 min.
Feline Corona Virus	Independent <sup>1</sup>	5 min.
Feline Immunodeficiency Virus	Independent <sup>1</sup>	5 min.
Hepatitis C	Independent <sup>1</sup>	5 min.
HIV	Independent <sup>1</sup>	5 min.
Norovirus (MVN)	EN 14476	1 min.
Parvovirus	EN 14476	10 min.
Poliovirus	EN 14476	1 min.
Rotavirus	EN 14476	1 min.
SARS	Independent <sup>1</sup>	5 min.

### Spores

<i>Bacillus subtilis</i>	EN 13704	30 min.
<i>Clostridium botulinum</i>	EN 1276	60 min.
<i>Clostridium difficile</i>	EN 13704	10 min.
<i>Clostridium perfringens</i>	EN 1276	15 min.

### Fungi

<i>Aspergillus niger</i>	EN 1275	5 min.
	EN 1650	5 min.
	EN 13697	15 min.
<i>Candida albicans</i>	EN 1275	1 min.
	EN 1650	5 min.
	EN 13624	1 min.
	EN 14562	15 min.

### Mycobacteria

<i>Mycobacterium avium</i>	EN 14348	5 min.
<i>Mycobacterium terrae</i>	EN 14348	5 min.

### Bacteria

<i>Campylobacter jejuni</i>	EN 1276	5 min.
	EN 13697	5 min.
<i>Corynebacterium ps.tuberculosis</i>	Independent <sup>2</sup>	10 min.
<i>Enterococcus hirae</i>	EN 1276	1 min.
	EN 13697	5 min.
	EN 13727	1 min.
	EN 14561	1 min.
<i>Escherichia coli</i>	EN 1276	1 min.
	EN 13697	5 min.
<i>Klebsiella pneumoniae</i>	EN 1276	1 min.
<i>Listeria monocytogenes</i>	EN 1276	1 min.
<i>Mannheimia haemolytica</i>	EN 1276	10 min.
Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA)	EN 1276	1 min.
<i>Proteus vulgaris</i>	EN 1276	5 min.
<i>Pseudomonas aeruginosa</i>	EN 1276	1 min.
	EN 13727	1 min.
	EN 14561	1 min.
<i>Salmonella typhimurium</i>	Independent <sup>2</sup>	10 min.
<i>Staphylococcus aureus</i>	EN 1276	1 min.
	EN 13697	5 min.
	EN 13727	1 min.
	EN 14561	1 min.
<i>Vibrio cholerae</i>	EN 1276	5 min.

### Other results

DEFRA approved under General Orders

Tested against bacteria over 7 days on stored solution

Independent testing by: <sup>1</sup> the University of Pretoria, RSA;

<sup>2</sup> Glasgow Caledonian University, Scotland, UK

Contact us for further details on our test results.